

# How to test for Fentanyl

## IV Users Guide

### Lower Your Risk of Overdose

#### Be Careful

- Go slow by taking smaller amounts
- Use with people you trust and inject at different times
- Have Naloxone available and know how to use it



### Step 1: Testing Cooker Residue

#### Prepare Shot


- Prepare your shot as normal and set aside the needle to test your substance left over in the cooker before injection
- Add enough water to fill your cooker. Stir if necessary



### Step 2: Testing

#### Use the Test Strip



- Place the test strip wavy side down in the water for 15s. Make sure the water is up to the line
- Remove strip from water and place it on a flat surface for 1-2 minutes



## Step 3: Results

### Read the Results

- A test **POSITIVE** for fentanyl is indicated by **ONE RED LINE**
- A test **NEGATIVE** for fentanyl is indicated by **TWO RED LINES**
- Remember, no test is 100% accurate, & a negative test does not necessarily mean you are safe



**NOTE:** It should be considered a **NEGATIVE** test even if a second red line is faded or faint  
If possible, use an additional test to confirm your results

